

SYNE2 Polyclonal Antibody

Cat No: HR1AP11000

For research use only

Overview

Product Name	SYNE2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	757kDa
GeneID?Human?	SYNE2 KIAA1011 NUA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	spectrin repeat containing nuclear envelope protein 2(SYNE2) Homo sapiens The protein encoded by this gene is a nuclear outer membrane protein that binds cytoplasmic F-actin. This binding tethers the nucleus to the cytoskeleton and aids in the maintenance of the structural integrity of the nucleus. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009],