

NOXO1 Polyclonal Antibody

Cat No: HR1AP10997

For research use only

Overview

Product Name	NOXO1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	NOXO1 P41NOX SH3PXD5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NADPH oxidase organizer 1(NOXO1) Homo sapiens This gene encodes an NADPH oxidase (NOX) organizer, which positively regulates NOX1 and NOX3. The protein contains a PX domain and two SH3 domains. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jun 2012],