

NDUA7 Polyclonal Antibody

Cat No: HR1AP10980

For research use only

Overview

Product Name	NDUA7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	NDUFA7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NADH:ubiquinone oxidoreductase subunit A7(NDUFA7) Homo sapiens This gene encodes a subunit of NADH:ubiquinone oxidoreductase (complex I), which is a multiprotein complex located in the inner mitochondrial membrane. Complex I functions in the transfer of electrons from NADH to the respiratory chain. [provided by RefSeq, Mar 2011],