

NDUA5 Polyclonal Antibody

Cat No: HR1AP10978

For research use only

Overview

Product Name	NDUA5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	NDUFA5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>NADH:ubiquinone oxidoreductase subunit A5(NDUFA5) Homo sapiens This nuclear gene encodes a conserved protein that comprises the B13 subunit of complex I of the mitochondrial respiratory chain. The encoded protein localizes to the inner mitochondrial membrane, where it is thought to aid in the transfer of electrons from NADH to ubiquinone. Alternative splicing results in multiple transcript variants. There are numerous pseudogenes of this gene on chromosomes 1, 3, 6, 8, 9, 11, 12, and 16. [provided by RefSeq, Apr 2014],</p>