

NDUAB Polyclonal Antibody

Cat No: HR1AP10975

For research use only

Overview

Product Name	NDUAB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	15kDa
GeneID?Human?	NDUFA11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NADH:ubiquinone oxidoreductase subunit A11(NDUFA11) Homo sapiens This gene encodes a subunit of the membrane-bound mitochondrial complex I. Complex I is composed of numerous subunits and functions as the NADH-ubiquinol reductase of the mitochondrial electron transport chain. Mutations in this gene are associated with severe mitochondrial complex I deficiency. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Oct 2010],