

NDUAA Polyclonal Antibody

Cat No: HR1AP10974

For research use only

Overview

Product Name	NDUAA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	NDUFA10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NADH:ubiquinone oxidoreductase subunit A10(NDUFA10) Homo sapiens The protein encoded by this gene is a component of 42 kDa complex I, the first enzyme complex in the electron transport chain of mitochondria. This protein has NADH dehydrogenase activity and oxidoreductase activity. It transfers electrons from NADH to the respiratory chain. A mutation in this gene was found in an individual with Leigh syndrome. [provided by RefSeq, Apr 2016],