

MY18B Polyclonal Antibody

Cat No: HR1AP10970

For research use only

Overview

Product Name	MY18B Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	282kDa
GeneID?Human?	MYO18B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>myosin XVIIIIB(MYO18B) Homo sapiens The protein encoded by this gene may regulate muscle-specific genes when in the nucleus and may influence intracellular trafficking when in the cytoplasm. The encoded protein functions as a homodimer and may interact with F actin. Mutations in this gene are associated with lung cancer. [provided by RefSeq, Jul 2008].</p>