

# MY18A Polyclonal Antibody

Cat No: HR1AP10969

For research use only

## Overview

Product Name	MY18A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	225kDa
GeneID?Human?	MYO18A KIAA0216 MYSPDZ
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:The TIAF1 protein is coded in the 3'-UTR region of MYO18A.,function:May be involved in the maintenance of the stromal cell architectures required for cell to cell contact.,similarity:Contains 1 IQ domain.,similarity:Contains 1 myosin head-like domain.,similarity:Contains 1 PDZ (DHR) domain.,subcellular location:Colocalizes with actin.,subcellular location:Lacks the PDZ domain. Diffusely localized in the cytoplasm.,subunit:Homodimer. Interacts with JAK3 and actin.,