

MYO5B Polyclonal Antibody

Cat No: HR1AP10964

For research use only

Overview

Product Name	MYO5B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	203kDa
GeneID?Human?	MYO5B KIAA1119
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	myosin VB(MYO5B) Homo sapiens The protein encoded by this gene, together with other proteins, may be involved in plasma membrane recycling. Mutations in this gene are associated with microvillous inclusion disease. [provided by RefSeq, Sep 2009],