

MYO3B Polyclonal Antibody

Cat No: HR1AP10963

For research use only

Overview

Product Name	MYO3B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	147kDa
GeneID?Human?	MYO3B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>myosin IIIB(MYO3B) Homo sapiens This gene encodes one of the class III myosins. Myosins are ATPases, activated by actin, that move along actin filaments in the cell. This class of myosins are characterized by an amino-terminal kinase domain and shown to be present in photoreceptors. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Mar 2014],</p>