

# MYO1B Polyclonal Antibody

Cat No: HR1AP10961

For research use only

## Overview

Product Name	MYO1B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	124kDa
GeneID?Human?	MYO1B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Motor protein that may participate in process critical to neuronal development and function such as cell migration, neurite outgrowth and vesicular transport.,similarity:Contains 1 myosin head-like domain.,similarity:Contains 6 IQ domains.,