

# MYBPH Polyclonal Antibody

Cat No: HR1AP10959

For research use only

## Overview

Product Name	MYBPH Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	MYBPH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function: Binds to myosin; probably involved in interaction with thick myofilaments in the A-band., similarity: Belongs to the immunoglobulin superfamily. MyBP family., similarity: Contains 2 fibronectin type-III domains., similarity: Contains 2 Ig-like C2-type (immunoglobulin-like) domains., tissue specificity: Skeletal muscle.,