

MYPC2 Polyclonal Antibody

Cat No: HR1AP10958

For research use only

Overview

Product Name	MYPC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	125kDa
GeneID?Human?	MYBPC2 MYBPCF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>myosin binding protein C, fast type(MYBPC2) Homo sapiens This gene encodes a member of the myosin-binding protein C family. This family includes the fast-, slow- and cardiac-type isoforms, each of which is a myosin-associated protein found in the cross-bridge-bearing zone (C region) of A bands in striated muscle. The protein encoded by this locus is referred to as the fast-type isoform. Mutations in the related but distinct genes encoding the slow-type and cardiac-type isoforms have been associated with distal arthrogyrosis, type 1 and hypertrophic cardiomyopathy, respectively. [provided by RefSeq, Jul 2012],</p>