

MYEOV Polyclonal Antibody

Cat No: HR1AP10952

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | MYEOV Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 34kDa |
| GeneID?Human? | MYEOV OCIM |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | developmental stage:Overexpressed in tumor cells lines with a t(11;14)(q13;q32) translocation., |