

## MLF2 Polyclonal Antibody

Cat No: HR1AP10950

For research use only

### Overview

Product Name	MLF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	MLF2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the MLF family.,tissue specificity:Ubiquitously expressed.,