

MEGF6 Polyclonal Antibody

Cat No: HR1AP10945

For research use only

Overview

Product Name	MEGF6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	169kDa
GeneID?Human?	MEGF6 EGFL3 KIAA0815
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 25 EGF-like domains.,