

MNX1 Polyclonal Antibody

Cat No: HR1AP10932

For research use only

Overview

Product Name	MNX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	MNX1 HLXB9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>motor neuron and pancreas homeobox 1(MNX1) Homo sapiens This gene encodes a nuclear protein, which contains a homeobox domain and is a transcription factor. Mutations in this gene result in Currarino syndrome, an autosomic dominant congenital malformation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009],</p>