

CDHR5 Polyclonal Antibody

Cat No: HR1AP10934

For research use only

Overview

Product Name	CDHR5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	CDHR5 MUCDHL MUPCDH UNQ2781/PRO7168
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cadherin related family member 5(CDHR5) Homo sapiens This gene is a novel mucin-like gene that is a member of the cadherin superfamily. While encoding nonpolymorphic tandem repeats rich in proline, serine and threonine similar to mucin proteins, the gene also contains sequence encoding calcium-binding motifs found in all cadherins. The role of the hybrid extracellular region and the specific function of this protein have not yet been determined. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jan 2010],</p>