

MOT11 Polyclonal Antibody

Cat No: HR1AP10926

For research use only

Overview

Product Name	MOT11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	SLC16A11 MCT11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:It is uncertain whether Met-1 or Met-25 is the initiator.,function:Proton-linked monocarboxylate transporter. Catalyzes the rapid transport across the plasma membrane of many monocarboxylates.,similarity:Belongs to the major facilitator superfamily. Monocarboxylate porter (TC 2.A.1.13) family.,