

# TAM41 Polyclonal Antibody

Cat No: HR1AP10924

For research use only

## Overview

Product Name	TAM41 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	TAMM41 C3orf31
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in the translocation of transit peptide-containing proteins across the mitochondrial inner membrane.,similarity:Belongs to the MMP37 family.,