

MOT1 Polyclonal Antibody

Cat No: HR1AP10925

For research use only

Overview

Product Name	MOT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	SLC16A1 MCT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 16 member 1(SLC16A1) Homo sapiens The protein encoded by this gene is a proton-linked monocarboxylate transporter that catalyzes the movement of many monocarboxylates, such as lactate and pyruvate, across the plasma membrane. Mutations in this gene are associated with erythrocyte lactate transporter defect. Alternatively spliced transcript variants have been found for this gene.[provided by RefSeq, Oct 2009],</p>