

MTFR1 Polyclonal Antibody

Cat No: HR1AP10920

For research use only

Overview

Product Name	MTFR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	MTFR1 CHPPR KIAA0009
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mitochondrial fission regulator 1(MTFR1) Homo sapiens This gene encodes a mitochondrial protein that is characterized by a poly-proline rich region. A chicken homolog of this protein promotes mitochondrial fission and the mouse homolog protects cells from oxidative stress. A related pseudogene of this gene is found on chromosome X. [provided by RefSeq, Mar 2009],