



# M1IP1 Polyclonal Antibody

Cat No: HR1AP10914

For research use only

## Overview

Product Name	M1IP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	M1D1IP1 MIG12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:It is uncertain whether Met-1 or Met-2 is the initiator.,function:Involved in stabilization of microtubules.,similarity:Belongs to the SPOT14 family.,subcellular location:Associated with microtubules.,