

# PGRC1 Polyclonal Antibody

Cat No: HR1AP10910

For research use only

## Overview

Product Name	PGRC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	PGRMC1 HPR6.6 PGRMC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	progesterone receptor membrane component 1(PGRMC1) Homo sapiens This gene encodes a putative membrane-associated progesterone steroid receptor. The protein is expressed predominantly in the liver and kidney. [provided by RefSeq, Mar 2010],