

MTSS1 Polyclonal Antibody

Cat No: HR1AP10911

For research use only

Overview

Product Name	MTSS1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	MTSS1 KIAA0429 MIM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>domain:The WH2 motif at the C-terminus binds to actin monomers.,function:May be related to cancer progression or tumor metastasis in a variety of organ sites, most likely through an interaction with the actin cytoskeleton.,similarity:Belongs to the MTSS1 family.,similarity:Contains 1 IMD (IRSp53/MIM homology) domain.,similarity:Contains 1 WH2 domain.,subunit:Binds to actin. Binds to the cytoplasmic domain of receptor protein tyrosine phosphatase delta.,tissue specificity:Expressed in many tissues, including spleen, thymus, prostate, testis, uterus, colon, and peripheral blood.,</p>