

# MRAP Polyclonal Antibody

Cat No: HR1AP10906

For research use only

## Overview

Product Name	MRAP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	MRAP C21orf61 FALP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	melanocortin 2 receptor accessory protein(MRAP) Homo sapiens This gene encodes a melanocortin receptor-interacting protein. The encoded protein regulates trafficking and function of the melanocortin 2 receptor in the adrenal gland. The encoded protein can also modulate signaling of other melanocortin receptors. Mutations in this gene have been associated with familial glucocorticoid deficiency type 2. Alternatively spliced transcript variants have been described. [provided by RefSeq, Dec 2009],