

MED4 Polyclonal Antibody

Cat No: HR1AP10904

For research use only

Overview

Product Name	MED4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	MED4 ARC36 DRIP36 VDRIP HSPC126
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mediator complex subunit 4(MED4) Homo sapiens This gene encodes a component of the Mediator complex. The Mediator complex interacts with DNA-binding gene-specific transcription factors to modulate transcription by RNA polymerase II. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012],