

# MED30 Polyclonal Antibody

Cat No: HR1AP10902

For research use only

## Overview

Product Name	MED30 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	MED30 THRAP6 TRAP25
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mediator complex subunit 30(MED30) Homo sapiens The multiprotein TRAP/Mediator complex facilitates gene expression through a wide variety of transcriptional activators. MED30 is a component of this complex that appears to be metazoan specific (Baek et al., 2002 [PubMed 11909976]).[supplied by OMIM, Nov 2010],