

MED29 Polyclonal Antibody

Cat No: HR1AP10901

For research use only

Overview

Product Name	MED29 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	MED29 IXL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mediator complex subunit 29(MED29) Homo sapiens MED29 is a subunit of the Mediator complex, a multiprotein coactivator of RNA transcription that interacts with DNA-bound transcriptional activators, RNA polymerase II (see MIM 180660), and general initiation factors (Sato et al., 2003 [PubMed 14576168]).[supplied by OMIM, Aug 2009],