

MED22 Polyclonal Antibody

Cat No: HR1AP10899

For research use only

Overview

Product Name	MED22 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	MED22 SURF5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mediator complex subunit 22(MED22) Homo sapiens This gene encodes a protein component of the mediator complex, which functions in the regulation of transcription by bridging interactions between gene-specific regulatory factors, RNA polymerase II, and general transcription factors. Alternatively spliced transcript variants encoding different isoforms have been observed. [provided by RefSeq, Jul 2013],