

MED10 Polyclonal Antibody

Cat No: HR1AP10896

For research use only

Overview

Product Name	MED10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	MED10 L6 TRG17 TRG20
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mediator complex subunit 10(MED10) Homo sapiens MED10 is a component of the Mediator complex, which is a coactivator for DNA-binding factors that activate transcription via RNA polymerase II (Sato et al., 2003 [PubMed 12584197]).[supplied by OMIM, Oct 2008],