

MATN2 Polyclonal Antibody

Cat No: HR1AP10895

For research use only

Overview

Product Name	MATN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	105kDa
GeneID?Human?	MATN2 UNQ193/PRO219
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>matrilin 2(MATN2) Homo sapiens This gene encodes a member of the von Willebrand factor A domain containing protein family. This family of proteins is thought to be involved in the formation of filamentous networks in the extracellular matrices of various tissues. This protein contains five von Willebrand factor A domains. The specific function of this gene has not yet been determined. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008].</p>