

MALD2 Polyclonal Antibody

Cat No: HR1AP10893

For research use only

Overview

Product Name	MALD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	MARVELD2 TRIC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>MARVEL domain containing 2(MARVELD2) Homo sapiens The protein encoded by this gene is a membrane protein found at the tight junctions between epithelial cells. The encoded protein helps establish epithelial barriers such as those in the organ of Corti, where these barriers are required for normal hearing. Defects in this gene are a cause of deafness autosomal recessive type 49 (DFNB49). Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011],</p>