



# LRP4 Polyclonal Antibody

Cat No: HR1AP10890

For research use only

## Overview

Product Name	LRP4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	209kDa
GeneID?Human?	LRP4 KIAA0816 LRP10 MEGF7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	LDL receptor related protein 4(LRP4) Homo sapiens This gene encodes a member of the low-density lipoprotein receptor-related protein family. The encoded protein may be a regulator of Wnt signaling. Mutations in this gene are associated with Cenani-Lenz syndrome. [provided by RefSeq, May 2010],