

# LRP1B Polyclonal Antibody

Cat No: HR1AP10888

For research use only

## Overview

Product Name	LRP1B Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	505kDa
GeneID?Human?	LRP1B LRPDIT
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	LDL receptor related protein 1B(LRP1B) Homo sapiens This gene encodes a member of the low density lipoprotein (LDL) receptor family. These receptors play a wide variety of roles in normal cell function and development due to their interactions with multiple ligands. Disruption of this gene has been reported in several types of cancer. [provided by RefSeq, Jun 2016],