

LRP12 Polyclonal Antibody

Cat No: HR1AP10887

For research use only

Overview

Product Name	LRP12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	94kDa
GeneID?Human?	LRP12 ST7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	LDL receptor related protein 12(LRP12) Homo sapiens This gene encodes a member of the low-density lipoprotein receptor related protein family. The product of this gene is a transmembrane protein that is differentially expressed in many cancer cells. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2010],