

LDH6B Polyclonal Antibody

Cat No: HR1AP10882

For research use only

Overview

Product Name	LDH6B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	LDHAL6B LDHAL6 LDHL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>catalytic activity:(S)-lactate + NAD(+) = pyruvate + NADH.,caution:It is uncertain whether Met-1 or Met-50 is the initiator.,developmental stage:Higher expression level in adult testis as compared to 6-week-old fetal testis.,pathway:Fermentation; pyruvate fermentation to lactate; (S)-lactate from pyruvate: step 1/1.,similarity:Belongs to the LDH/MDH superfamily. LDH family.,tissue specificity:Testis specific.,</p>