

LDHC Polyclonal Antibody

Cat No: HR1AP10883

For research use only

Overview

Product Name	LDHC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	LDHC LDH3 LDHX
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	lactate dehydrogenase C(LDHC) Homo sapiens Lactate dehydrogenase C catalyzes the conversion of L-lactate and NAD to pyruvate and NADH in the final step of anaerobic glycolysis. LDHC is testis-specific and belongs to the lactate dehydrogenase family. Two transcript variants have been detected which differ in the 5' untranslated region. [provided by RefSeq, Jul 2008],