

# LDH6A Polyclonal Antibody

Cat No: HR1AP10881

For research use only

## Overview

Product Name	LDH6A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	LDHAL6A LDHL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:(S)-lactate + NAD(+) = pyruvate + NADH.,pathway:Fermentation; pyruvate fermentation to lactate; (S)-lactate from pyruvate: step 1/1.,similarity:Belongs to the LDH/MDH superfamily.,similarity:Belongs to the LDH/MDH superfamily. LDH family.,