

LIN7A Polyclonal Antibody

Cat No: HR1AP10878

For research use only

Overview

Product Name	LIN7A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	LIN7A MALS1 VELI1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	lin-7 homolog A, crumbs cell polarity complex component(LIN7A) Homo sapiens The protein encoded by this gene is involved in generating and maintaining the asymmetric distribution of channels and receptors at the cell membrane. The encoded protein also is required for the localization of some specific channels and can be part of a protein complex that couples synaptic vesicle exocytosis to cell adhesion in the brain. [provided by RefSeq, May 2016],