

LRBA Polyclonal Antibody

Cat No: HR1AP10879

For research use only

Overview

Product Name	LRBA Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	314kDa
GeneID?Human?	LRBA BGL CDC4L LBA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>LPS responsive beige-like anchor protein(LRBA) Homo sapiens The protein encoded by this gene is a member of the WDL-BEACH-WD (WBW) gene family. Its expression is induced in B cells and macrophages by bacterial lipopolysaccharides (LPS). The encoded protein associates with protein kinase A and may be involved in leading intracellular vesicles to activated receptor complexes, which aids in the secretion and/or membrane deposition of immune effector molecules. Defects in this gene are associated with the disorder common variable immunodeficiency-8 with autoimmunity. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2012].</p>