

LHX4 Polyclonal Antibody

Cat No: HR1AP10875

For research use only

Overview

Product Name	LHX4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	LHX4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	LIM homeobox 4(LHX4) Homo sapiens This gene encodes a member of a large protein family which contains the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein is a transcription factor involved in the control of differentiation and development of the pituitary gland. Mutations in this gene cause combined pituitary hormone deficiency 4. [provided by RefSeq, Dec 2010],