

LRC42 Polyclonal Antibody

Cat No: HR1AP10874

For research use only

Overview

Product Name	LRC42 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	LRRC42
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the LRRC42 family.,similarity:Contains 2 LRR (leucine-rich) repeats.,