

LARP1 Polyclonal Antibody

Cat No: HR1AP10873

For research use only

Overview

Product Name	LARP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	120kDa
GeneID?Human?	LARP1 KIAA0731 LARP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the LARP family.,similarity:Contains 1 HTH La-type RNA-binding domain.,