

## A4GAT Polyclonal Antibody

Cat No: HR1AP10872

For research use only

### Overview

Product Name	A4GAT Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	A4GALT A14GALT A4GALT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>alpha 1,4-galactosyltransferase(A4GALT) Homo sapiens The protein encoded by this gene catalyzes the transfer of galactose to lactosylceramide to form globotriaosylceramide, which has been identified as the P(k) antigen of the P blood group system. This protein, a type II membrane protein found in the Golgi, is also required for the synthesis of the bacterial verotoxins receptor. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Dec 2015].</p>