

UBP54 Polyclonal Antibody

Cat No: HR1AP10863

For research use only

Overview

Product Name	UBP54 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	185kDa
GeneID?Human?	USP54 C10orf29
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:Although the active site residues are conserved, it lacks the conserved His residue which is normally found 9 residues before the catalytic His.,function:Has no peptidase activity.,similarity:Belongs to the peptidase C19 family.,tissue specificity:Weakly expressed in a few tissues.,