

HOME1 Polyclonal Antibody

Cat No: HR1AP10856

For research use only

Overview

Product Name	HOME1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	HOMER1 SYN47
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	homer scaffolding protein 1(HOMER1) Homo sapiens This gene encodes a member of the homer family of dendritic proteins. Members of this family regulate group 1 metabotropic glutamate receptor function. [provided by RefSeq, Jul 2008],