

TF2LY Polyclonal Antibody

Cat No: HR1AP10855

For research use only

Overview

Product Name	TF2LY Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	TGIF2LY TGIFLY
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>TGFB induced factor homeobox 2 like, Y-linked(TGIF2LY) Homo sapiens This gene encodes a member of the TALE/TGIF homeobox family of transcription factors. This gene lies within the male specific region of chromosome Y, in a block of sequence that is thought to be the result of a large X-to-Y transposition. The C-terminus of this protein is divergent from that of its chromosome X homolog (TGIF2LX), suggesting that this protein may act as a regulator of TGIF2LX. [provided by RefSeq, Jul 2008],</p>