



SIX4 Polyclonal Antibody

Cat No: HR1AP10852

For research use only

Overview

Product Name	SIX4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	85kDa
GeneID?Human?	SIX4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SIX homeobox 4(SIX4) Homo sapiens This gene encodes a member of the homeobox family, subfamily SIX. The drosophila homolog is a nuclear homeoprotein required for eye development. Studies in mouse show that this gene product functions as a transcription factor, and may have a role in the differentiation or maturation of neuronal cells. [provided by RefSeq, May 2010],