

NKX62 Polyclonal Antibody

Cat No: HR1AP10849

For research use only

Overview

Product Name	NKX62 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	NKX6-2 GTX NKX6B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 homeobox DNA-binding domain.,tissue specificity:Highest expression in brain.,